

Biotin Mouse IgG2a, κ Isotype Control[C1.18.4]

Catalog No. E-AB-F09803B
Storage Store at 2~8°C, Avoid freeze / thaw cycles

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Details

Form	Liquid
Concentration	0.5 mg/mL
Size	25µg/100µg
Clone No.	C1.18.4
Host	Mouse
Isotype	Mouse IgG2a, κ
Application	FCM
Storage Buffer	Phosphate buffered solution, pH 7.2, containing 0.09% sodium azide and 1% BSA.
Shipping	Biological ice pack at 4 °C
Stability & Storage	Keep as concentrated solution. Store at 2~8°C .Do not freeze. This product is guaranteed up to one year from purchase.

Recommended usage

Each lot of this antibody is quality control tested by flow cytometric analysis as negative control. Use at concentrations comparable to those of the specific antibody of interest.

Related Information

1. 流式细胞术样本制备技术 <https://www.elabscience.cn/List-detail-10425.html>
2. 细胞内细胞因子染色步骤 <https://www.elabscience.cn/List-detail-10424.html>
3. 流式实验常见问题及解决方案 <https://www.elabscience.cn/List-detail-5690.html>
4. 如何通过光谱图选择合适的检测通道? <https://www.elabscience.cn/List-detail-16891.html>