

Ki67 Polyclonal Antibody

Catalog Number:E-AB-40528



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Zebrafish
Immunogen	Recombinant protein corresponding to Zebrafish Ki67.
Host	Rabbit
Isotype	IgG
Purification	Antigen Affinity Purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol,pH 7.4

Applications Recommended Dilution

IF	1:100-200
-----------	-----------

Data



Immunofluorescence of Zebrafish eyeball with anti Ki67 Polyclonal Antibody at dilution of 1:100.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Ki67 is a nuclear protein that is expressed in proliferating cells. Ki67 is preferentially expressed during late G1, S, M, and G2 phases of the cell cycle, while cells in the G0 (quiescent) phase are negative for this protein. Increased proliferativ

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Tel: 400-999-2100

Email: techsupport@elabscience.cn

Web: www.elabscience.cn