

ANXA5 Polyclonal Antibody

Catalog Number:E-AB-40114



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant human Annexin A5 protein
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.05% Proclin300, 50% glycerol, pH7.3.

Applications Recommended Dilution


WB	1:500-1:1000
IHC	1:50-1:100

Data




Western Blot analysis of HeLa and Mouse kidney using ANXA5 Polyclonal Antibody at dilution of 1:600

Observed Mw:36KDa
Calculated Mw:36KDa



Immunohistochemistry of paraffin-embedded Rat kidney using ANXA5 Polyclonal Antibody at dilution of 1:50



Immunohistochemistry analysis of paraffin-embedded Human Prostate using ANXA5 Polyclonal Antibody(Elabscience® Product Detected by Lifespan).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Tel: 400-999-2100

Email: techsupport@elabscience.cn

Web: www.elabscience.cn

ANXA5 Polyclonal Antibody

Catalog Number:E-AB-40114



involved in the blood coagulation cascade.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Tel: 400-999-2100

Email: techsupport@elabscience.cn

Web: www.elabscience.cn