

NOS3 Monoclonal Antibody

Catalog Number:E-AB-22088



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG
Clone	Clone:6G3
Purification	Protein A purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.4

Applications Recommended Dilution

WB 1:500-2000

Data



www.elabscience.cn

Western Blot analysis of Rat heart using NOS3
Monoclonal Antibody at dilution of 1:1000.

Observed Mw:130kDa
Calculated Mw:133kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Produces nitric oxide (NO) which is implicated in vascular smooth muscle relaxation through a cGMP-mediated signal transduction pathway. NO mediates vascular endothelial growth factor (VEGF)-induced angiogenesis in coronary vessels and promotes blood clotting through the activation of platelets.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Tel: 400-999-2100

Email: techsupport@elabscience.cn

Web: www.elabscience.cn