

# Rabbit Anti-Goat IgG(H+L)(peroxidase/HRP conjugated)



Catalog Number: E-AB-1004 3 Publications

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Immunogen</b>	Goat IgG
<b>Host</b>	Rabbit
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	HRP
<b>Formulation</b>	0.01M PBS with 50% glycerol, PH7.4

## Applications Recommended Dilution

<b>Western Blotting</b>	1:8000-15000(Enhanced chemiluminescent detection)
<b>Western Blotting</b>	1:1000-5000(DAB detection)
<b>Direct ELISA</b>	1:4000-8000(TMB detection)
<b>Immunohisto/cytochemistry</b>	1:400-4000(DAB detection)

## Preparation & Storage

**Storage** Store at -20°C. Avoid repeated freeze / thaw cycles.

## For Research Use Only