

# Goat Anti-Rabbit IgG(H+L)(Cyanine3 conjugated)

Catalog Number:E-AB-1010 7 Publications



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Immunogen</b>	Rabbit IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	CY3 (Ex/Em:554nm/568-574nm)
<b>Formulation</b>	0.01M PBS with 50% glycerol and 0.05% Sodium Azide PH7.6

## Applications Recommended Dilution

<b>Immunofluorescence</b>	1:50-100
---------------------------	----------

## Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid repeated freeze / thaw cycles.Protected from prolonged exposure to light.
----------------	---

**For Research Use Only**

Tel: 400-999-2100

Web: [www.elabscience.cn](http://www.elabscience.cn)

Email: [techsupport@elabscience.cn](mailto:techsupport@elabscience.cn)