

Goat Anti-Rabbit IgG(H+L)(FITC conjugated)

Catalog Number:E-AB-1014 17 Publications



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Immunogen	Rabbit IgG
Host	Goat
Purification	Affinity purification
Conjugation	FITC (Ex/Em:495nm/525nm)
Formulation	0.01M PBS with 50% glycerol and 0.05% Sodium Azide PH7.6

Applications Recommended Dilution

Immunofluorescence	1:20-100
---------------------------	----------

Preparation & Storage

Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.Protected from prolonged exposure to light.
----------------	---

For Research Use Only

Tel: 400-999-2100

Web: www.elabscience.cn

Email: techsupport@elabscience.cn