

# Goat Anti-Mouse IgG(H+L)(Elab Fluor<sup>®</sup> 488 conjugated)



Catalog Number: E-AB-1056 6 Publications

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Immunogen</b>	Mouse IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Elab Fluor <sup>®</sup> 488 (Ex/Em:493nm/519nm)
<b>Formulation</b>	0.01M PBS with 50% glycerol and 0.05% Sodium Azide PH7.6

## Applications Recommended Dilution

<b>Immunofluorescence</b>	1:100-300
---------------------------	-----------

## Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid repeated freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	--

**For Research Use Only**

Tel: 400-999-2100

Web: [www.elabscience.cn](http://www.elabscience.cn)

Email: [techsupport@elabscience.cn](mailto:techsupport@elabscience.cn)