

Goat Anti-Mouse IgG(H+L)(peroxidase/HRP conjugated)

Catalog Number:E-AB-1001 55 Publications

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Immunogen	Mouse IgG
Host	Goat
Purification	Affinity purification
Conjugation	HRP
Formulation	0.01M PBS with 50% glycerol, PH7.4

Applications Recommended Dilution

Western Blotting	1:5000-50000 (Enhanced chemiluminescent detection)
Western Blotting	1:1000-5000(DAB detection)
Direct ELISA	1:5000-30000(TMB detection)
Immunohisto/cytochemistry,1:400-4000(DAB detection)	

Preparation & Storage

Storage	Store at -20°C. Avoid repeated freeze / thaw cycles.
----------------	--

For Research Use Only