

# Goat Anti-Mouse IgG(H+L)(peroxidase/HRP conjugated)



Catalog Number: E-AB-1001 55 Publications

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Immunogen</b>	Mouse IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	HRP
<b>Formulation</b>	0.01M PBS with 50% glycerol, PH7.4

## Applications Recommended Dilution

<b>Western Blotting</b>	1:5000-50000 (Enhanced chemiluminescent detection)
<b>Western Blotting</b>	1:1000-5000(DAB detection)
<b>Direct ELISA</b>	1:5000-30000(TMB detection)
<b>Immunohisto/cytochemistry,1:400-4000(DAB detection)</b>	

## Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid repeated freeze / thaw cycles.
----------------	--

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Tel: 400-999-2100

Email: [techsupport@elabscience.cn](mailto:techsupport@elabscience.cn)

Web: [www.elabscience.cn](http://www.elabscience.cn)