

mCherry Monoclonal Antibody

Catalog Number:E-AB-20087

2 Publications



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG
Clone	Clone:2F4
Purification	Protein A purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.4

Applications Recommended Dilution

WB 1:5000-1:10000

Data



Western Blot analysis of mCherry recombinant protein using mCherry Monoclonal Antibody at dilution of 1) 1:5000 2) 1:10000.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Tel: 400-999-2100

Email: techsupport@elabscience.cn

Web: www.elabscience.cn